EN

Konelab™ / T Series HbA1c Pretreatment Liquid

REF 981924 12 x 20 ml

**THIS PACKAGE INSERT IS APPLICABLE FOR USE
OUTSIDE THE US. ANY REFERENCE TO THE KONELAB
SYSTEMS ALSO REFERS TO THE T SERIES.**

INTENDED USE

For on-board and manual hemolyzing of whole blood samples for *in vitro* determination of hemoglobin A1c (HbA1c) on Konelab analyzers. Refer to 981658 package insert for details on the HbA1c test. All test results must be interpreted with regard to the clinical context.

PRINCIPLE OF THE PROCEDURE (1)

This method uses TTAB (tetradecyl trimethyl ammonium bromide) as the detergent in the hemolyzing reagent to eliminate interference from leucocytes (TTAB does not lyse leucocytes). Sample pretreatment to remove labile HbA1c is not necessary.

REAGENT INFORMATION

HbA1c Pretreatment Liquid 12 x 20 ml

Concentrations

Tetradecyl trimethyl ammonium bromide (TTAB) 9 g/l

Precautions

For *in vitro* diagnostic use only. Exercise the normal precautions required for handling all laboratory reagents.

Preparation

Reagent is ready for use.

Storage and Stability

Reagents in unopened vials are stable at 2...8 °C until the expiry date printed on the label. Refer to the Application Notes of your Konelab analyzer for the on-board stability of reagents. Once opened the reagents are stable for 30 days if stored closed at 2...8 °C when not in use and if contamination is avoided.

SPECIMEN COLLECTION AND PREPARATION**Sample Type**

Capillary, EDTA- or heparinized blood. (1,2)
The sample is stable for (3) 3 days at 15...25 °C
7 days at 2...8 °C
6 months at -20 °C
(freeze once only).

Precautions

Human samples should be handled and disposed of as if they were potentially infectious.

Sample preparation and stability**1) On-board hemolyzing**

HbA1c Pretreatment liquid is inserted to the reagent register of the analyzer.

Note 1: Check that there are no bubbles in the bottleneck or on the surface of the reagent when you insert the reagent vials or vessels in the Konelab analyzer.

Note 2: The reagent solution must be in the reagent disk temperature of analyzer before performing the on board hemolyzing.

The samples are inserted to the instrument in the violet coloured pretreatment segment (codes 984800, 984801) in primary tubes. Immediately prior to inserting the sample to the instrument, mix the sample thoroughly to obtain uniform distribution of the erythrocytes. Avoid the formation of foam. Start analysing immediately after inserting the segment to the instrument. Insert up to 28 samples (two full segments) at a time. If using Konelab 20, insert a maximum of 14 samples (one full segment) at a time. Further samples can be inserted after the results of the first 14 are ready.

The whole blood samples should be run as a separate batch to avoid possible carry over and minimize sample stand-in-time in the instrument.

Samples that appear to have a high erythrocyte sedimentation rate should be put in the first positions (1-4) of the pretreatment segment. Samples with very high sedimentation rate should be hemolyzed manually. Also samples with abnormally high viscosity (very high hemoglobin or hematocrit) should be hemolyzed manually.

During the on-board hemolyzing the sample is diluted 1+100 with the HbA1c Pretreatment Liquid. If another dilution ratio is required, the hemolyzing must be carried out manually.

A maximum of 5 HbA1c requests can be given for one sample. When multiple requests are given the instrument does one hemolyzate from which the requests are made.

When using on-board hemolyzing it is important to carry out the daily maintenance procedures regularly (cleaning the dispensers, mixers, washing wells and dispensing table).

NOTE! Only whole blood samples in primary tubes are inserted in a violet pretreatment segment to the instrument. **Do not insert whole blood samples in a normal black segment.** Controls, calibrators and manually hemolyzed samples are inserted in a normal black segment to the instrument. **Do not insert controls, calibrators or manually hemolyzed samples in a violet pretreatment segment.**

2) Manual hemolyzing

(Refer to 981658 HbA1c package insert)

TEST PROCEDURE

Refer to the Reference Manual and Application Notes for an on-board hemolyzing procedure on your Konelab analyzer. Any application which has not been validated by Thermo Fisher Scientific Oy cannot be performance guaranteed and therefore must be evaluated by the user.

Refer to 981658 HbA1c package insert for details.

Materials provided

HbA1c Pretreatment Liquid, code 981924, 12 x 20 ml

Materials required but not provided

HbA1c reagent kit , code 981658, calibrators included.
HbA1c Control Normal, code 981695 and HbA1c Control Abnormal, code 981696
0.9 % NaCl
Pretreatment segment, codes 984800, 984801.

BIBLIOGRAPHY

- Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers, Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

MANUFACTURED BY

Thermo Fisher Scientific Oy
Clinical Diagnostics Finland
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland
Tel. +358 9 329 100, Fax +358 9 3291 0300
www.thermo.com/konelab

Date of revision (yyyy-mm-dd)

2009-03-03

Changes from previous version

Pretreatment segment code 984801 added.

DE

Konelab™ / T Series HbA1c Pretreatment Liquid

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Tetradecyltrimethylammoniumbromid (TTAB) 9 g/l

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- Bakkere DL et al. Multicenter Evaluation of an Improved Immunturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers. Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results. 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical variables. Broschüre in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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BIBLIOGRAPHIE

- Bakkere DL et al. Multicenter Evaluation of an Improved Immunturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers. Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results. 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical Variables. Broschüre in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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Koncentrace

Tetradecyltrimethylamonium-bromid (TTAB) 9 g/l

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- 1 Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers. Clin Lab 1999;45:123-137.
- 2 Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- 3 Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GfT Verlag 1996.

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ET
Konelab™ / T Series
HbA1c Pretreatment Liquid

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KIRJANDUS

1. Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers, Clin Lab 1999;45:123-137.
2. Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
3. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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- Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers. Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results, 1ª edición, TH-Books Verlagsgesellschaft, Fráncfort/Main, Alemania, 1998: 142-148.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Folleto incluido en: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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Thermo Fisher Scientific Oy
Clinical Diagnostics Finland
Ratastie 2, P.O. Box 100, FI-01621 Vantaa Finland
Tel. +358 9 329 100, Fax +358 9 3291 0300
www.thermo.com/konelab

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BIBLIOGRÁFIA

- Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers. Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition. TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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IT
Konelab™ / T Series
HbA1c Pretreatment Liquid

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BIBLIOGRAFIA

- Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers. Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition. TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- Guder WG, Narayanan S, Wissner H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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LITERATŪRA

1. Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers, Clin Lab 1999;45:123-137.
2. Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
3. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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1. Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers, Clin Lab 1999;45:123-137.
2. Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory Results, 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main, Germany, 1998, 142-148.
3. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996

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1. Bakkeren DL et al. Multicenter Evaluation of an Improved Immunturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers, Clin Lab 1999;45:123-137.
2. Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
3. Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GIT Verlag 1996.

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BIBLIOGRAFIA

- Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers, Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition, TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. Darmstadt: GIT Verlag 1996.

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Thermo Fisher Scientific Oy
Clinical Diagnostics Finland
Ratastie 2, P.O. Box 100, FI-01621 Vantaa, Finland
Tel.: +358 9 329 100, Fax +358 9 3291 0300
www.thermo.com/konelab

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BIBLIOGRAFIE

- Bakkeren DL et al. Multicenter Evaluation of an Improved Immunoturbidimetric Assay for the Determination of HbA1c on Clinical Chemistry Analyzers. Clin Lab 1999;45:123-137.
- Thomas L (ed.) Clinical Laboratory Diagnostics: Use and Assessment of Clinical Laboratory results, 1st edition. TH-Books Verlagsgesellschaft, Frankfurt/Main Germany 1998: 142-148.
- Guder WG, Narayanan S, Wisser H, Zawta B. List of Analytes; Preanalytical Variables. Brochure in: Samples: From the patient to the laboratory. Darmstadt: GFT Verlag 1996.

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