

Eliminate Specimen Identification Errors in Cytology Slide Processing

Thermo Scientific SlideMate Slide Printer

Simply scan the bar code on the specimen vial and SlideMate™ prints a slide with the correct information for that specimen. This eliminates the opportunity for error when matching slides and vials at the point of slide processing. Each slide is printed individually, on demand, so that only one slide—the right one—is available to load with the vial.

- On-demand slide printing reduces errors and improves productivity
- Prints directly onto the slide eliminating the need for labels or handwriting
- Small footprint fits easily into the Cytology preparation area
- Chemical resistant media ensures print reliability

To learn more about on-demand slide printing at the point of Cytology slide preparation, call:

800-522-7270 ext.564

or visit: www.thermo.com/slidemate



"Slide labeling alone accounted for 67% of defects"

Source: The Henry Ford Production Systems Defects and Waste in Surgical Pathology as a Basis for Quality Improvement Initiatives.

©2007 American Society for Clinical Pathology

On-demand slide printing at the point of Cytology preparation is just another way we support you—the face behind the glass.

SlideMate Benefits for your Laboratory

- On-demand slide printing streamlines workflow and increases accuracy of identification
- Small footprint fits easily into the lab
- Flexible data input methods make for easy integration into the laboratory
- High-quality printing resistant to harsh chemicals
- Slides are held inside the unit's side door to protect them from dust and the environment
- User can select the number of slides to print per case, which enhances productivity in the lab

Specifications

Width	6.0", 152mm
Height	8.9", 226mm
Depth	10.4", 264mm
Weight	4.3kg, 9.5lbs
Power Requirements	Universal 100-240 VAC 47-63 Hz 1.2A Maximum
Color	Sterile Steel
Slide Type	3" x 1" or 25mm x 75mm x 1mm with frosted or colored end
Slide Capacity	150 pieces (input stack) 10 pieces (output tray)
Media Type	Thermal transfer ribbon, 1" wide roll
Media Capacity	1 roll, approximately 12,000 Slides
Media Color	Black
Throughput	6 Slides/minute typical



Order Information

Description	Order Number
SlideMate Print On-Demand Printing System Includes SlideMate Printer, Scanner, Keyboard Interface, Mini Keyboard, Box of SlidesD	B8130001
SlideMate with LabWriter™ Software Includes SlideMate Printer, LabWriter Basic Software for SlideMate, Box of Slides	B8130002
SlideMate with PC – Premium Package Includes SlideMate Printer, All-in-one PC, Digital Imager, LabWriter Basic Software, Box of Slides	B8130003
SlideMate Printer Includes SlideMate Printer	B8130004
PrintMate Available Accessories	Order Number
Digital Imager	B35310006
All-in-one PC	A79510021
LabWriter Basic Software for SlideMate	B81310001
Network Cable	B35310012
Crossover Ethernet Cable	B81310005
Keyboard Interface	B81320043
Mini Keyboard	B81310004
Thermal Transfer Ribbon	B81320041

© 2009 Thermo Fisher Scientific Inc. All rights reserved. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Anatomical Pathology

4481 Campus Drive
Kalamazoo, MI 49008
(800) 522-7270

Robert-Bosch-Str. 49
D-69190 Walldorf/Germany
+ 49(0)6227-8360

93/96 Chadwick Road
Astmoor, Runcorn Cheshire WA7 1 PR
800 018 9396, + 1 (0) 1928 562501